June 15, 2004

For Research Use Only. Not for use in diagnostic procedures.

Detection of caspase activity

APOPCYTO Caspase-3 Colorimetric Assay Kit

CODE No. 4800
100 assays

Before use, thoroughly read these Instructions.

Description

Caspase is a member of the cysteine aspartic acid-specific protease family, which is activated by a variety of signals of death receptor ligation, DNA damages, serum starvation and stresses *etc.* Active caspase recognizes a lot of several molecules as substrates to cleave them, occurring to biological events corresponding to the apoptosis. For example, ICAD (inhibitor of caspase-activated deoxyribonuclease) is inactivated and CAD (caspase-activated deoxyribonuclease) is indirectly activated by caspase-3, and it is related to chromatin fragmentation for nucleosome units. Caspase recognizes several structural proteins as a substrate to cleave them, and the cleavage is associated with the unique apoptosis cell morphology of chromatin condensation, nucleus fragmentation and cytoplasmic integrity. Active caspase could act in apoptosis process, so that detection of caspase activity is necessary for apoptosis research frequently. Though there are numerous ways to measure of caspase activity, colorimetric or fluorometric substrates are used popularly.

Assay Principle

Active caspase recognizes 4 amino acid residues in a substrate molecule, which are from the aspartate residue next to N-terminal 3 amino acid residues, and specifically cleaves at the C-terminal side of the aspartate residue. So this characterization is useful for detecting caspase activity by using synthetic substrates. In this method, 4 amino acid sequences is labeled with *p*NA (*p*-nitroanilide), MCA (4-Methyl Coumaryl 7-amide) or AFC (7-amino-4-trifluoromethyl coumarin) at the C-terminal side. Free *p*NA, AMC(7-amino-4-methyl-coumarin) or AFC is released from the labeled synthetic substrate on cleavage by active caspase, and monitored by a 96 well microplate reader (*p*NA: absorbance wavelength 400 or 405 nm, AMC: excitation wavelength 380 nm, fluorescence wavelength 460 nm, AFC: excitation wavelength 400 nm, fluorescence wavelength 505 nm). For example, if caspase is active in sample, tetrapeptide-*p*NA as a synthetic substrate is cleaved and the wavelength is shifted as shown in Fig.1. Because the amount of free *p*NA is proportional to the amount of caspase activity present in the sample, caspase activity can be calculated by monitoring of the optical density of free *p*NA at wavelength 405nm.

APOPCYTO Caspase-3 Colorimetric Assay Kit is provided to detect caspase-3 activity in cell extract with DEVD-*p*NA as a substrate, which DEVD sequence is recognized by active caspase-3 selectively. APOPCYTO Caspase-3 Colorimetric Assay Kit is composed of reagents necessary for the detection, which are Cell Lysis buffer, Reaction buffer and so on. Caspase-3 activity can be detected easily and quickly using APOPCYTO Caspase-3 Colorimetric Assay Kit.

CODE No.4800



Fig. 1 Comparison of UV-Absorption Spectrum of Peptidyl-pNA with that of pNA

Intended Use

For research use only. Not for use in *in vitro* diagnostic procedures for clinical diagnosis.

Materials	Quantity	Color
Cell Lysis Buffer	50 mL X 1 vial	Clear
2X Reaction Buffer	2 mL X 3 vials	Blue
1M DTT	500 μ L X 1 vial	Yellow
Substrate DEVD-p NA (10 mM)	550 μ L X 1 vial	Brown
Inhibitor DEVD-FMK (1 mM)	55 μ L X 1 vial	Orange
100 mM p NA (p -nitroanilide)	200 μ L X 1 vial	Green

Product Components

Expiration

Please see the label of this kit.

Storage Conditions

Store at -20°C. After thaw, the Cell Lysis Buffer and 2X Reaction Buffer can be stored at 4°C.

Material to Be Supplied by the User

- microcentrifuge (For harvest cells)
- 96-well microplate reader
(For measure the absorbance in the well at 400-405 nm)
- micropipette
- 96-well microplate (flat-bottom, clear polystyrene)
- microcentrifuge tube

Note

- When caspase-3 activity is not measured immediately after cell preparations, store the cell lysates or cell pellets at -20°C and use them within 1 week.
- Store caspase-3 substrate away from light.
- Ensure that DTT is added to 2X Reaction Buffer to obtain a final concentration of 10 mM before use.

Preparation of reagents

- 1) Add DTT to 2X Reaction Buffer to obtain final concentration of 10 mM just before use. (For example, 10 µL of 1 M DTT to 990 µL of 2X Reaction Buffer)
- 2) Preparation of *p*NA (*p*-nitroanilide) standards
 - i) Dilute 100 mM *p*NA (*p*-nitroanilide) to 5 mM *p*NA in Cell Lysis Buffer (5 µL of 100 mM *p*NA is added to 95 µL of Cell Lysis Buffer).
 - ii) Dilute 5 mM *p*NA in Cell Lysis Buffer according to the following table to give these final concentrations.

<i>p</i>NA concentration	<i>p</i>NA	Cell Lysis Buffer
500 µM (50 nmole)	50 µL of 5 mM <i>p</i> NA	450 µL
250 µM (25 nmole)	250 µL of 500 µM <i>p</i> NA	250 µL
125 µM (12.5 nmole)	250 µL of 250 µM <i>p</i> NA	250 µL
62.5 µM (6.25 nmole)	250 µL of 125 µM <i>p</i> NA	250 µL
31.25 µM (3.125 nmole)	250 µL of 62.5 µM <i>p</i> NA	250 µL
15.625 µM (1.5625 nmole)	250 µL of 31.25 µM <i>p</i> NA	250 µL
0 µM (0 nmole)		250 µL

Caspase Assay Protocol

1. Induce apoptosis in cells using the desired method, Concurrently incubate a control culture without induction.
2. Count the cells and harvest the $1-5 \times 10^6$ cells by centrifugation at 400 X *g* for 5 minutes.
3. Resuspend the cells in ice-cold Cell Lysis Buffer (50-500 µL), and

incubate the cells on ice for 10 minutes.

4. Centrifuge the cell lysates at 10,000 X *g* for 5 minutes at 4°C to precipitate cellular debris. Transfer the supernatants (cell extracts) to new microcentrifuge tubes and put on ice.
If necessary, total protein concentration of the cell extracts may be measured by standard protein assay methods, and adjust protein concentration to 100-200 µg / 50 µL in Cell Lysis Buffer.
5. Add 50 µL of 2X Reaction Buffer containing 10 mM DTT to each well of a 96 well microplate.
10 mM DTT should be added to 2X Reaction Buffer just before use.
6. Add 50 µL of cell lysates to each well. Prepare wells which were added 50 µL of Cell Lysis Buffer in stead of cell lysates to measure blank absorbance.
7. [Optional method]
To verify that the signal detected by the kit is due to protease activity, incubate an induced sample with DEVD-FMK before adding substrate.
Add 50 µL of cell extracts to 50 µL of 2X Reaction Buffer (containing 10 mM DTT) and 1 µL of 1 mM DEVD-FMK, and sequentially add 5 µL of Caspase-3 Substrate. Then, incubate at 37°C for 1-2 hours, and measure the absorbance with other samples.
8. Add 5 µL of Caspase-3 Substrate to each well.
9. Cover the plate with a plate sealer and incubate at 37 °C for 1-2 hours (or up to overnight maximum).
10. Add 100 µL of *p*NA (*p*-nitroanilide) standards to empty wells.
11. Measure the absorbance in the wells at 400 or 405 nm. Calculate the specific activity (SA) of caspase-3 present in each sample using following formula:
 - i) Construct a standard curve using the absorbance of *p*NA (*p*-nitroanilide) standards.
 - ii) Calculate the relative absorbance (Ar) as:
 $A_a = (\text{induced apoptosis sample } A_{405}) - (\text{blank } A_{405})$
 $A_n = (\text{negative sample } A_{405}) - (\text{blank } A_{405})$
 $A_r = A_a - A_n$
 - iii) Calculate liberated free *p*NA concentration (B µM) equal to Ar by *p*NA standard curve.
 - iv) Calculate caspase-3 activity (C) present in each sample as:
C = nmole *p*NA liberated per hour
= $B \times 0.1 \text{ mL} / \text{incubation hour}$
= $0.1B / \text{incubation hour}$
 - v) If total protein concentration of the cell extracts is D mg/mL, Calculate the specific activity (SA) of caspase-3 present in each sample as:
SA = C per µg protein
= $(0.1B / \text{incubation hour}) / (D \text{ mg/mL} \times 0.1 \text{ mL})$
= $B / (D \cdot \text{incubation hour})$

Assay Example


Fig. 2 Caspase-3 activity in CH-11 treated Jurkat cell lysates After Jurkat cells were incubated with 100 ng/mL of anti-Fas monoclonal antibody (CH-11) in the presence or absence of caspase-3 inhibitor DEVD-FMK at 37 °C for 4 hours, caspase-3 activity was measured by APOPCYTO Caspase-3 Colorimetric Assay Kit.

References

- 1) Wolf, B. B., *et al.*, *J. Biol. Chem.* **274**, 20049-20052 (1999)
- 2) Thornberry, N. A., *et al.*, *J. Biol. Chem.* **272**, 17907-17911 (1997)
- 3) Walker, N. P., *et al.*, *Cell* **78**, 343-352 (1994)
- 4) Sleath, P. R., *et al.*, *J. Biol. Chem.* **265**, 14526-14528 (1990)

Manufacture

MBL MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.

URL <https://ruo.mbl.co.jp>

e-mail support@mbi.co.jp, TEL 052-238-1904

Before use, thoroughly read these Instructions.

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

Fig. 1 Comparison of UV-Absorption Spectrum of Peptidyl-*p*NA with that of *p*NA

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