

Hypoxia Probe / Hypoxia Probe Solution

Product Name	Quantity	Cat. #	MW	Storage
Hypoxia Probe	2 mg	LOX-1 NC-LOX-1P (For Asia)	711.87	Room temp., Protection from light
Hypoxia Probe Solution	2 mM, 100 μ L	LOX-1S NC-LOX-1S (For Asia)		-20°C, Protection from light

Please refer to expiration date on the label.


Precautions for use

Thank you for purchasing a Hypoxia probe. Please read all the instructions before use.

Introduction

This probe is a phosphorescent light-emitting iridium complex. Phosphorescence of this probe is quenched by oxygen, and it is increased in response to low levels of oxygen which is detectable by a general fluorescent microscopy (red-fluorescence). The signal change is reversible, and this probe easily permeates cell membrane. So, the reagent shows real-time response according to oxygen levels in living cells. This probe can be used to observe hypoxia condition inside spheroid and detect subcutaneous xenograft tumors in small animals.

Spectral Characteristics



Absorption spectrum is given in green and phosphorescence in red. Absorption and phosphorescence are peaked at 483 and 616 nm, respectively. ¹⁾

Example of use for Hypoxia Probe / Hypoxia Probe Solution

1. Hypoxia Probe LOX-1 is dissolved in 1.4 mL DMSO to make 2 mmol / L stock solution. (LOX-1 stock solution is stored at -20°C except for use.)

Note: Please use fresh or good-condition DMSO that not absorbed moisture.

2. The stock solution of Hypoxia Probe or Hypoxia Probe solution is diluted with culture



medium to prepare 20 $\mu\text{mol/L}$ working solution just before use (100 times dilution).

3. Then, Add 10 μL -aliquots of the working solution into 100 μL culture medium gently (final concentration of the probe is 2 $\mu\text{mol/L}$).
4. After incubation for one day, red-phosphorescence is observed with a general fluorescent microscopy with the standard TRITC filter set (for example G-2A filter block: Ex 510-560, DM575, BA590).

Note: Incubation time is depend on your sample. For monolayer cells, incubation time may be enough 0.5 hrs. For spheroids, it may be 8-12 hrs.

Note: We recommend using FITC excitation and Texas Red[®] emission filters for best results.

References

1. Zhang S, Hosaka M, Yoshihara T, Negishi K, Iida Y, Tobita S, Takeuchi T. Phosphorescent Light-Emitting Iridium Complexes Serve as a Hypoxia-Sensing Probe for Tumor Imaging in Living Animals. *Cancer Res.* 2010 Jun 1;70(11):4490-8.

Product warranty

All materials supplied in product pass an inspection at ORGANOGENIX, Inc., prior to shipment. However, if you received any defective product, please contact us.

CONTACT INFORMATION

Sales and Customer Support:

Japan/Asia: MEDICAL & BIOLOGICAL LABORATORIES CO., LTD. (MBL)

Email: support@mbi.co.jp

URL: ruo.mbl.co.jp, (Japanese), URL: ruo.mbl.co.jp/bio/e/index.html, (English)

Product URL, Japanese: www.nanocultureplate.jp, English: www.nanocultureplate.com

Tel: 052-238-1904

China: MBL Beijing Biotech Co., Ltd.

Email: info@bio-med.com.cn

URL: www.bio-med.com.cn/Home.php

Tel: +86-10-82899503

Americas, Europe and other regions: MBL International Corporation

Email: sales@mblintl.com

URL: www.mblintl.com/

Tel: +1-800-200-5459, Fax: +1-781-939-6963