

# Anti-Hemoglobin F (Human) pAb

<b>CODE No.</b>	PM078
<b>CLONALITY</b>	Polyclonal
<b>ISOTYPE</b>	Rabbit Ig, affinity purified
<b>QUANTITY</b>	100 µL
<b>SOURCE</b>	Purified IgG from rabbit serum
<b>IMMUNOGEN</b>	KLH conjugated synthetic peptide, CLTSLGDAIKH (corresponding to amino acid residues 69-78 of human hemoglobin, gamma A)
<b>REACTIVITY</b>	This antibody reacts with hemoglobin F, and does not cross-react with hemoglobin A.
<b>FORMULATION</b>	PBS containing 1% BSA and 0.09% NaN <sub>3</sub>
	*Azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush plenty of water when disposing materials containing azide into drain.
<b>STORAGE</b>	This antibody solution is stable for one year from the date of purchase when stored at 4°C.

## APPLICATION-CONFIRMED

Immunohistochemistry 1:100 (paraffin section)

\*Heat or enzymatic treatment is required.

## SPECIES CROSS REACTIVITY on IHC

Species	Human	Mouse	Rat	Hamster
Tissues	Fetal liver, MDS bone marrow	Not tested	Not tested	Not tested
Reactivity	+			

**Entrez Gene ID** 3047 (Human)

**REFERENCE** 1) Choi, J. W., *et al.*, *Int. J. Hematol.* **74**, 277-280 (2001)

For more information, please visit our website at <https://ruo.mbl.co.jp/>.

The descriptions of the following protocols are examples. Each user should determine the appropriate condition.


### **Immunohistochemistry for formalin fixed paraffin-embedded section**

*Immunohistochemistry was performed using Leica BOND-MAX™ and Bond™ reagents. Leica Bond-MAX™ and BOND™ reagents are trademarks of Leica Biosystems Newcastle Ltd. For more information, please contact Leica Biosystems Newcastle Ltd.*


- 1) Soak slides in Deparaffinized Solution for 30 min. at 72°C.
- 2) Wash the slides 3 times with Ethanol.
- 3) Wash the slides 1 time in Washing Solution for 5 min.
- 4) Perform heat or enzymatic treatment for antigen retrieval.  
Heat-treatment: Soak slides in Epitope Retrieval Solution for 20 min. at 100°C and then let the slides cool down to room temperature for 12 min. at room temperature.  
Enzymatic treatment: Soak slides in Enzymatic Solution for 10 min. at 37°C.
- 5) Wash the slides 3 times in Washing Solution.
- 6) Incubate the slides with primary antibody diluted with Antibody Diluent as suggested in the **APPLICATION** for 15 min. at room temperature. (The concentration of antibody will depend on the conditions.)
- 7) Wash the slides 4 times in Washing Solution.
- 8) Incubate with Polymer Enhancer reagent for 8 min. at room temperature.
- 9) Wash the slides 4 times in Washing Solution.
- 10) Incubate with Polymer reagent for 10 min. at room temperature.
- 11) Wash the slides 2 times in Washing Solution.
- 12) Wash the slides 2 times in distilled water.
- 13) Immerse the slides in Blocking reagent for 5 min. at room temperature.
- 14) Wash the slides 4 times in distilled water.
- 15) Visualize by reacting for 10 min. with DAB reagent. \*DAB is a suspect carcinogen and must be handled with care. Always wear gloves.
- 16) Wash the slides 4 times in distilled water.
- 17) Counterstain with Hematoxylin for 3 min., wash the slides in distilled water.
- 18) Dehydrate by immersing in Ethanol, followed by immersing in Deparaffinized Solution. Now ready for mounting.

(Positive controls for Immunohistochemistry; Bone marrow from patient with MDS and autopsied fetal liver.)


**MDS bone marrow**



**Adult bone marrow**



**Fetal liver**



### ***Immunohistochemical detection of Hemoglobin F***

Brown: Anti-Hemoglobin F (Human) pAb (MBL, code no. PM078)  
Blue: Hematoxylin

Data were kindly provided by Dr. Masafumi Ito.  
(Japanese Red Cross Nagoya Daiichi Hospital)