

**RiboCluster Profiler™**

RBP Antibody

# Anti-TARDBP (TDP-43) pAb

<b>CODE No.</b>	RN107PW
<b>CLONALITY</b>	Polyclonal
<b>ISOTYPE</b>	Rabbit Ig, affinity purified
<b>QUANTITY</b>	100 µL, 1 mg/mL
<b>SOURCE</b>	Purified Ig from rabbit serum
<b>FORMURATION</b>	PBS containing 50% Glycerol (pH 7.2). No preservative is contained.
<b>STORAGE</b>	This antibody solution is stable for one year from the date of purchase when stored at -20°C.

## APPLICATIONS

<u>Western blotting</u>	1:500 for chemiluminescence detection system
<u>Immunoprecipitation</u>	5 µL/500 µL of nuclear extract from 1 x 10 <sup>7</sup> cells/sample

## SPECIES CROSS REACTIVITY on WB

Species	Human	Mouse	Rat	Hamster
Cell	HeLa, HEK293T, Jurkat, K562	NIH/3T3, WR19L	Rat1	CHO
Reactivity	+	+	-	+

**Entrez Gene ID** 23435 (Human), 230908 (Mouse)

For more information, please visit our web site <http://ruo.mbl.co.jp/je/rip-assay/>

**LICENSING OPPORTUNITY:** The RIP-Assay uses patented technology (US patent No. 6,635,422, US patent No. 7,504,210, JP patent No. 5,002,105) of Ribonomics, Inc. MBL manufactures and distributes this product under license from Ribonomics, Inc. Researchers may use this product for their own research. Researchers are not allowed to use this product or RIP-Assay technology for commercial purpose without a license. For commercial use, please contact us for licensing opportunities at [RIP@mbi.co.jp](mailto:RIP@mbi.co.jp)

## RELATED PRODUCTS

### RIP-Assay Kit

RN1001	RIP-Assay Kit
RN1005	RIP-Assay Kit for <i>microRNA</i>

### RIP-Certified Antibody

RN001P	Anti-EIF4E pAb
RN002P	Anti-EIF4G1 pAb
RN003P	Anti-EIF4G2 pAb
RN004P	Anti-ELAVL1 (HuR) pAb
RN005P	Anti-ELAVL2 (HuB) pAb
RN006P	Anti-ELAVL3 (HuC) pAb
RN007P	Anti-IGF2BP1 (IMP1) pAb
RN008P	Anti-IGF2BP2 (IMP2) pAb
RN009P	Anti-IGF2BP3 (IMP3) pAb
RN010P	Anti-MSI1 (Musashi1) pAb
RN011P	Anti-PTBP1 pAb
RN012P	Anti-STAU1 pAb
RN013P	Anti-STAU2 pAb
RN014P	Anti-TIA1 pAb
RN015P	Anti-YBX1 pAb
RN016P	Anti-FMR1 pAb
RN017P	Anti-FXR1 pAb
RN018P	Anti-FXR2 pAb
RN019P	Anti-HNRNPK pAb
RN020P	Anti-ILF3 pAb
RN021P	Anti-KHDRBS1 pAb
RN022P	Anti-PABPC4 pAb
RN024P	Anti-PCBP1 pAb
RN025P	Anti-PCBP2 pAb
RN026P	Anti-PUM1 pAb
RN027P	Anti-PUM2 pAb
RN028P	Anti-EIF2C1 (AGO1) pAb
RN032P	Anti-CIRBP pAb
RN033P	Anti-TNRC6A (GW182) pAb
RN037P	Anti-AUH pAb
RN038P	Anti-CPEB1 pAb
RN041P	Anti-KHDRBS2 (SLM1) pAb
RN045P	Anti-SLBP pAb
RN001M	Anti-IGF2BP1 (IMP1) mAb (6H6)
RN003M	Anti-EIF2C2 (AGO2) mAb (1B1-E2H5)
RN004M	Anti-Ribosomal P0/P1/P2 mAb (9D5)
RN005M	Anti-EIF2C2 (AGO2) mAb (2A8)
RN006M	Anti-EIF4E mAb (C107-3-5)
RN007M	Anti-ELAVL1 (HuR) mAb (C67-1)
RN009M	Anti-PABPC1 mAb (10E10)

### RBP Antibody

RBP Antibody works on WB and /or IP, but not certified for working on RIP-Assay.

RN008MW	Anti-ELAVL1 (HuR) mAb (C54-6)
RN023PW	Anti-PABPN1 pAb
RN031PW	Anti-ZFP36 (TTP) pAb
RN046PW	Anti-SYNERIP (HNRNPQ) pAb
RN047PW	Anti-PTBP2 pAb
RN048PW	Anti-G3BP1 pAb
RN049PW	Anti-G3BP2 pAb
RN050PW	Anti-GRSF1 pAb
RN051PW	Anti-HDLBP (Vigilin) pAb

RN052PW	Anti-HNRNPC pAb
RN053PW	Anti-PAIP1 pAb
RN054PW	Anti-PCBP3 pAb
RN056PW	Anti-SERBP1 pAb
RN057PW	Anti-TARBP1 pAb
RN058PW	Anti-TARBP2 pAb
RN059PW	Anti-TIAL1 pAb
RN060PW	Anti-HNRNPD (AUF1) pAb
RN061PW	Anti-HNRNPA0 pAb
RN062PW	Anti-DGCR8 pAb
RN063PW	Anti-DHX9 pAb
RN064PW	Anti-FUSIP1 (SRSF10) pAb
RN065PW	Anti-KHSRP pAb
RN066PW	Anti-KIAA0020 pAb
RN067PW	Anti-PPP1R10 pAb
RN068PW	Anti-PPP1R8 pAb
RN069PW	Anti-RBM14 pAb
RN077PW	Anti-SMN1 pAb
RN078PW	Anti-SMNDC1 pAb
RN079PW	Anti-SRSF7 (9G8) pAb
RN080PW	Anti-SRSF3 (SRp20) pAb
RN081PW	Anti-SRSF9 (SRp30c) pAb
RN082PW	Anti-SRSF5 (SRP40) pAb
RN084PW	Anti-SRRM1 (SRM160) pAb
RN085PW	Anti-U2AF1 pAb
RN086PW	Anti-U2AF2 pAb
RN087PW	Anti-ALYREF (THOC4) pAb
RN088PW	Anti-NXF1 (TAP) pAb
RN089PW	Anti-MAGOH pAb
RN090PW	Anti-DDX21 pAb
RN091PW	Anti-DDX23 pAb
RN092PW	Anti-NONO (P54NRB) pAb
RN093PW	Anti-PRPF4 pAb
RN094PW	Anti-PRPF8 pAb
RN095PW	Anti-SNRNP200 pAb
RN096PW	Anti-SNRNP40 pAb
RN097PW	Anti-SNRNP70 pAb
RN098PW	Anti-EDC4 pAb
RN099PW	Anti-EIF4A1 pAb
RN100PW	Anti-EXOSC5 (RRP46) (Human) pAb
RN101PW	Anti-FBL (Fibrillarlin) pAb
RN102PW	Anti-GEMIN2 (Human) pAb
RN103PW	Anti-NCBP1 (CBP80) pAb
RN104PW	Anti-PAN2 (USP52) (Human) pAb
RN105PW	Anti-PARN pAb
RN106PW	Anti-SFPQ (PSF) pAb
RN107PW	Anti-TARDBP (TDP-43) pAb
RN108PW	Anti-UPF1 pAb
RN109PW	Anti-XRN1 (Human) pAb


D216-3	Anti-hnRNP-A2/B1 mAb (C20308)
M162-3	Anti-p62 (SQSTM1) (Human) mAb (5F2)
PM045	Anti-p62 (SQSTM1) pAb
PM066	Anti-p62 C-terminal pAb

For the latest information of RiboCluster Profiler™, please visit our website at <http://ruo.mbl.co.jp/je/rip-assay/>

### **SDS-PAGE & Western blotting**

- 1) Wash  $1 \times 10^7$  cells 3 times with PBS and suspend them in 1 mL of Laemmli's sample buffer, then sonicate briefly (up to 20 sec.).
- 2) Boil the samples for 3 min. and centrifuge. Load 10  $\mu$ L of the sample per lane in a 1-mm-thick SDS-polyacrylamide gel (12.5% acrylamide) for electrophoresis.
- 3) Blot the protein to a polyvinylidene difluoride (PVDF) membrane at 1 mA/cm<sup>2</sup> for 1 hr. in a semi-dry transfer system (Transfer Buffer: 25 mM Tris, 190 mM glycine, 20% MeOH). See the manufacturer's manual for precise transfer procedure.
- 4) To reduce nonspecific binding, soak the membrane in 5% skimmed milk (in PBS, pH 7.2) for overnight at 4°C
- 5) Wash the membrane with PBS-T (0.05% Tween-20 in PBS) [5 min. x 3 times].
- 6) Incubate the membrane with primary antibody diluted with 1% skimmed milk (in PBS, pH 7.2) as suggested in the **APPLICATIONS** for 1 hr. at room temperature. (The concentration of antibody will depend on the conditions.)
- 7) Wash the membrane with PBS-T (10 min. x 3 times).
- 8) Incubate the membrane with the 1:5,000 of anti-IgG (Rabbit) pAb-HRP (MBL; code no. 458) diluted with 1% skimmed milk (in PBS, pH 7.2) for 1 hr. at room temperature.
- 9) Wash the membrane with PBS-T (10 min. x 3 times).
- 10) Wipe excess buffer on the membrane, then incubate it with appropriate chemiluminescence reagent for 1 min. Remove extra reagent from the membrane by dabbing with paper towel, and seal it in plastic wrap.
- 11) Expose to an X-ray film in a dark room for 3 min. Develop the film as usual settings. The condition for exposure and development may vary.

(Positive controls for Western blotting; HeLa, HEK93T, Jurkat, K562, NIH/3T3, WR19L and CHO)



#### ***Western blot analysis of TARDBP (TDP-43)***


- Lane 1: HeLa
- Lane 2: HEK293T
- Lane 3: Jurkat
- Lane 4: K562
- Lane 5: NIH/3T3
- Lane 6: WR19L
- Lane 7: Rat1
- Lane 8: CHO

Immunoblotted with Anti-TARDBP (TDP-43) pAb (RN107PW)

## **Immunoprecipitation**

- 1) Wash  $2 \times 10^7$  cells 4 times with PBS and resuspend them with 1 mL of ice-cold Lysis Buffer (+) (MBL; code no. RN1001) containing appropriate protease inhibitors and DTT. Vortex thoroughly, then incubate on ice for 5 min.
- 2) Centrifuge the tube at 12,000 x g for 5 min. at 4°C and discard the supernatant.
- 3) Wash the pellet 3 times with PBS and resuspend them with 500 µL of RIPA buffer, then sonicate briefly.
- 4) Centrifuge the tube at 12,000 x g for 5 min. at 4°C. and transfer the supernatant to another tube.
- 5) Add 500 µL of ice-cold Lysis buffer (+) into the supernatant. Mix well by pipetting up and down.
- 6) Add 20 µL of 50% protein G agarose beads slurry resuspended in ice-cold Lysis Buffer (+) into the supernatant. Incubate it at 4°C with rotating for 1 hr.
- 7) Centrifuge the tube at 2,000 x g for 2 min. at 4°C and transfer the supernatant to another tube (precleared sample).
- 8) Mix 20 µL of 50% protein G agarose beads slurry resuspended in 1 mL of ice-cold Wash Buffer (+) (MBL; code no. RN1001) containing DTT at the appropriate concentration with normal rabbit IgG (RIP-Assay Kit) or anti-TARDBP (TDP-43) pAb (MBL; code no. RN107PW) as suggested in the **APPLICATIONS**. Incubate at 4°C with rotating for 1 hr.
- 9) Wash the beads 1 time with ice-cold Lysis Buffer (+). Carefully discard the supernatant.
- 10) Add 500 µL of the precleared sample (prepared in step 7)) to the tube containing antibody conjugated beads, then incubate with gentle agitation for 2 hr. at 4°C.
- 11) Wash the beads 4 times with 1 mL of ice-cold Wash Buffer (+).
- 12) Resuspend the beads in 20 µL of Laemmli's sample buffer, boil for 3 min. and centrifuge.
- 13) Load 10 µL of the sample per lane in a 1-mm-thick SDS-polyacrylamide gel (12.5% acrylamide) for electrophoresis.
- 14) Blot the protein to a polyvinylidene difluoride (PVDF) membrane at 1 mA/cm<sup>2</sup> for 1 hr. in a semi-dry transfer system (Transfer Buffer: 25 mM Tris, 190 mM glycine, 20% MeOH). See the manufacturer's manual for precise transfer procedure.
- 15) To reduce nonspecific binding, soak the membrane in 5% skimmed milk (in PBS, pH 7.2) for overnight at 4°C.
- 16) Wash the membrane with PBS-T (0.05% Tween-20 in PBS) [5 min. x 3 times].
- 17) Incubate the membrane with primary antibody diluted with 1% skimmed milk (in PBS, pH 7.2) as suggested in the **APPLICATIONS** for 1 hr. at room temperature. (The concentration of antibody will depend on the conditions.)
- 18) Wash the membrane with PBS-T (10 min. x 3 times).
- 19) Incubate the membrane with 1: 1,000 of Rabbit TrueBlot<sup>®</sup> anti-Rabbit IgG-HRP (eBioscience; code no. 18-8816-33) diluted with 1% skimmed milk (in PBS, pH 7.2) for 1 hr. at room temperature.
- 20) Wash the membrane with PBS-T (10 min. x 3 times).
- 21) Wipe excess buffer on the membrane, then incubate it with appropriate chemiluminescence reagent for 1 min.
- 22) Remove extra reagent from the membrane by dabbing with paper towel, and seal it in plastic wrap.
- 23) Expose to an X-ray film in a dark room for 3 min. Develop the film as usual. The condition for exposure and development may vary.

(Positive control for Immunoprecipitation; HeLa nuclear extract)



### ***Immunoprecipitation of TARDBP (TDP-43) from HeLa nuclear extract***

Lane 1: Input  
Lane 2: Normal Rabbit IgG  
Lane 3: Anti-TARDBP (TDP-43) pAb (RN107PW)

Immunoblotted with RN107PW