

## Protocol for Fluorescent *in situ* hybridization (FISH) on cultured cells

### Day 1

#### *Cell seeding*

- ↓ Seed the proper number of cells on the dish (Nunc Lab-Tek II Chamber Slide, code no. 154534; Thermo Scientific) and incubate it at 37°C overnight.

### Day 2

#### *Pretreatment*

- ↓ Wash wells with PBS twice.
- ↓ Fix in 4% paraformaldehyde at room temperature for 20 min.
- ↓ Wash wells with PBS twice.
- ↓ Incubate cells in 0.2 N HCl at room temperature for 20 min.
- ↓ Wash wells with Nuclease free water.
- ↓ Incubate cells in 37°C prewarmed Proteinase K solution (0.6 µg/mL) at 37°C for 7 min. (Optimizing the Proteinase K concentration and the reaction time is recommended.)
- ↓ To stop the protease reaction, incubate cells in 0.2% Glycine/ PBS at room temperature for 10 min.
- ↓ Wash wells with PBS twice.
- ↓ Wash wells with Nuclease free water.
- ↓ Fix in 4% paraformaldehyde at room temperature for 20 min.
- ↓ Wash wells with PBS.
- ↓ Wash wells with Nuclease free water.
- ↓ Incubate cells in Acetylation solution (1.5% Triethanol amine, 0.03 N HCl, 0.25% Acetic anhydride) at room temperature for 15 min.
- ↓ Wash wells with PBS.
- ↓ Incubate cells in 55°C prewarmed Prehybridization solution (50% formamide, 1×Denhardt's solution, 2×SSC, 10 mM EDTA (pH8.0), 100 µg/mL yeast tRNA, 0.01% Tween-20) at 55°C for 2 hr.

#### *Hybridization*

- ↓ Denature 50 µg/sample of Probe solution (1 µg/mL DIG-labeled RNA probe, 400 µg/mL yeast tRNA diluted with formamide) at 80°C for 10 min., then quench at 4°C for 5 min.
- ↓ Mix 50 µL of 55°C prewarmed Hybridization Solution (2×Denhardt's solution, 4×SSC, 20 mM EDTA (pH8.0), 0.02% Tween-20, 10% Dextran sulfate) and 50 µL

of Probe Solution, then vortex well.

- ↓ Add the mixed solution and transfer the dish in humidity chamber to avoid dry. Incubate it at 55°C for at least 16 hr. (The condition of the hybridization reaction depends on the probe being used.)

### **Day 3**

#### *Wash step after hybridization*

- ↓ Wash wells twice with 55°C prewarmed Wash Buffer A (50% formamide, 2×SSC, 0.01% Tween 20) at 55°C for 30 min.
- ↓ Wash wells with NTET Buffer (10 mM Tris-HCl (pH8.0), 1 mM EDTA, 0.01% Tween 20, 500 mM NaCl).
- ↓ Incubate cells in 37°C prewarmed NTET Buffer with RNase A (10 µg/mL) at 37°C for 1 hr.
- ↓ Wash wells with NTET Buffer, twice.
- ↓ Wash wells with 55°C prewarmed Wash Buffer B (2×SSC, 0.01% Tween 20) at 55°C for 30 min.
- ↓ Wash wells twice with 55°C prewarmed Wash Buffer C (0.1×SSC, 0.01% Tween 20) at 55°C for 30 min.
- ↓ Wash wells with PBS.

#### *Detection*

- ↓ Add Blocking Buffer (1× Blocking Reagent (code no. 1158576200; Roche), 1×PBS, 0.01% Tween 20), and incubate cells at room temperature for 1 hr.
- ↓ Add Anti-Digoxigenin (DIG) mAb (MBL; code no. M227-3) diluted with Blocking Buffer at room temperature for 1 hr.
- ↓ Wash wells three times with 0.2% Tween 20/ PBS at room temperature for 5 min. each.
- ↓ Wash wells with PBS.
- ↓ Add Alexa fluor 594 conjugated anti-mouse IgG (code no. A11032; Invitrogen) diluted with Blocking Buffer at room temperature for 1 hr.
- ↓ Wash wells three times with 0.2% Tween 20/ PBS at room temperature for 5 min. each.
- ↓ Wash wells with PBS.
- ↓ Mount the slide with Mount medium with DAPI.