

**PEPTIDE**

# HA-tag peptide

Code No.	Quantity	Form
3320-205	2 mg x 5 vials	Lyophilized

**DESCRIPTION:**

This synthetic peptide including HA-tag sequence has been designed for the elution of HA-tagged Protein Purification Gel. The peptide sequence is YPYDVDPDYA.

**MOLECULAR WEIGHT:**

1,102 Da

**SUPPLIED:**

Lyophilized in PBS

Reconstitute in 1.0 mL of distilled water. After reconstitution, the HA-tag peptide concentration will be 2 mg/mL in PBS.

**STORAGE:**

Lyophilized HA-tag peptide is stable for one year from the date of shipment when stored at 4°C. After reconstitution, prepare appropriate aliquots and freeze them at -20°C. Avoid repeated freeze/thaw cycles.

**RELATED PRODUCTS:**

Please visit our website at <http://ruo.mbl.co.jp/>