

MONOCLONAL ANTIBODY

Anti-CX₃CR1 (Human) mAb-FITC

Code No.	Clone	Subclass	Quantity	Concentration
D070-4	2A9-1	Rat IgG2b κ	100 μL	500 μg/mL

BACKGROUND: There are several subfamilies in the chemokine superfamily. In addition to CXC, CC, and C subfamily, Fractalkine (FKN), which has the novel CX₃C chemokine motif and the mucin-like domain, has recently been identified and reported. This mucin-chemokine hybrid type of protein can exist in two forms; either membrane-bound form or soluble secreted form. The membrane-bound form of FKN protein is markedly induced on primary endothelial cells by inflammatory cytokines, and it promotes strong adhesion of NK cells and CD8⁺ T cells. The soluble secreted form of FKN can be released, presumably by proteolysis at a membrane-proximal dibasic cleavage site, and has chemotactic activity for these leukocytes.

CX₃CR1, which is recently identified FKN receptor, is also G-protein-coupled seven-transmembrane receptor as another chemokine receptor families, and is expressed on the cell surface of NK cells and CD8⁺ T cells. It is also reported that CX₃CR1 and FKN mediate both leukocytes migration and adhesion.

SOURCE: This antibody was purified from hybridoma (clone 2A9-1) supernatant using protein G agarose. This hybridoma was established by fusion of mouse myeloma cell P3U1 with WKY/NCrj rat lymph nodes immunized with non-mammalian cells expressing human CX₃CR1 protein.

FORMULATION: 50 μg IgG in 100 μL volume of PBS containing 1% BSA and 0.09% NaN₃.

*Azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush plenty of water when disposing materials containing azide into drain.

STORAGE: This antibody solution is stable for one year from the date of purchase when stored at 4°C.

REACTIVITY: This antibody reacts with human CX₃CR1 on Flow cytometry.

APPLICATION:

Flow cytometry; 10 μg/mL (final concentration)
*Please refer to the data sheet (MBL code no. D070-3) for other applications.

Detailed procedure is provided in the following **PROTOCOLS.**

SPECIES CROSS REACTIVITY:

Species	Human	Mouse	Rat
Cells	Transfectant, peripheral blood lymphocytes, monocyte	Not tested	Not tested
Reactivity on FCM	+		


INTENDED USE:

For Research Use Only. Not for use in diagnostic procedures.

REFERENCES:

- 1) Tae Yu, H., *et al.*, *Cell Mol. Immunol.* **12**, 466-473 (2015)
- 2) Liu, Y., *et al.*, *Circ. J.* **74**, 2585-2591 (2010)
- 3) Momose, K., *et al.* *Int. J. Hematol.* **86**, 174-179 (2007)
- 4) Bjerkeli, V., *et al.*, *Rheumatology* **46**, 1422-1427 (2007)
- 5) EI-Shazly, A., *et al.* *J. Immunol.* **176**, 1860-1868 (2006)
- 6) Lee, R. H., *et al.*, *Blood* **107**, 2153-2161 (2006)
- 7) Makishima H *et al.* *Leukemia* **19**, 1169-1174 (2005)
- 8) Vitale, M., *et al.*, *Eur. J. Immunol.* **34**, 1715-1722 (2004)

As this antibody is really famous all over the world, a lot of researches have been reported. These references are a part of such reports.



Flow cytometric analysis of CX₃CR1 expression on transfectant. Open histogram indicates the reaction of isotopic control to the cells. Shaded histogram indicates the reaction of D070-4 to the cells.

PROTOCOLS:

Flow cytometric analysis for floating cells


We usually use Fisher tubes or equivalents as reaction tubes for all steps described below.

- 1) Wash the cells 3 times with washing buffer [PBS containing 2% fetal calf serum (FCS) and 0.09% NaN₃].
- 2) Resuspend the cells with washing buffer (5x10⁶)

cells/mL).

- 3) Add 50 μ L of the cell suspension into each tube, and centrifuge at 500 x g for 1 minute at room temperature (20~25°C). Remove supernatant by careful aspiration.
- 4) Add 20 μ L of Clear Back (human Fc receptor blocking reagent, MBL; code no. MTG-001) to the cell pellet after tapping. Mix well and incubate for 5 minutes at room temperature.
- 5) Add 20 μ L of the primary antibody at the concentration as suggested in the **APPLICATIONS** diluted with the washing buffer. Mix well and incubate for 30 minutes at room temperature.
- 6) Add 1 mL of the washing buffer followed by centrifugation at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 7) Resuspend the cells with 500 μ L of the washing buffer and analyze by a flow cytometer.

(Positive control for Flow cytometry; Transfectant)



Flow cytometric analysis of CX₃CR1 expression on peripheral blood lymphocyte and monocyte. The staining intensity of D070-4 is shown in the horizontal axis with CD16 staining on the vertical axis.

Flow cytometric analysis for whole blood cells

We usually use Fisher tubes or equivalents as reaction tubes for all steps described below.

- 1) Add 50 μ L of the primary antibody at the concentration as suggested in the **APPLICATIONS** diluted with the washing buffer into each tube.
- 2) Add 50 μ L of whole blood into each tube. Mix well, and incubate for 30 minutes at room temperature (20~25°C).

- 3) Add 1 mL of washing buffer [PBS containing 2% fetal calf serum (FCS) and 0.09% NaN₃] followed by centrifugation at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 4) Lyse with OptiLyse C (for analysis on Beckman Coulter instruments) or OptiLyse B (for analysis on BD instruments), using the procedure recommended in the respective package inserts.
- 5) Add 1 mL of H₂O to each tube and incubate for 10 minutes at room temperature.
- 6) Centrifuge at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 7) Add 1 mL of washing buffer followed by centrifugation at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 8) Resuspend the cells with 500 μ L of the washing buffer and analyze by a flow cytometer.

(Positive controls for Flow cytometry; Human lymphocyte and monocyte)

RELATED PRODUCTS:

- D070-3 Anti-CX₃CR1 (Human) mAb (2A9-1)
- D070-5 Anti-CX₃CR1 (Human) mAb-PE (2A9-1)
- D070-A48 Anti-CX₃CR1 (Human) mAb-Alexa Fluor[®] 488 (2A9-1)
- D063-3 Anti-CD191 (CCR1) (Human) mAb (#141-2)
- D063-5 Anti-CD191 (CCR1) (Human) mAb-PE (#141-2)
- D085-3 Anti-CD193 (CCR3) (Human) mAb (444-11)
- D085-4 Anti-CD193 (CCR3) (Human) mAb-FITC (444-11)
- D085-5 Anti-CD193 (CCR3) (Human) mAb-PE (444-11)
- D074-3 Anti-CD197 (CCR7) (Human) mAb (6B3)
- D124-3 Anti-CD195 (CCR5) (Human) mAb (T227)
- D124-4 Anti-CD195 (CCR5) (Human) mAb-FITC (T227)
- D123-3 Anti-CD184 (CXCR4) mAb (A145)
- D123-4 Anti-CD184 (CXCR4) mAb-FITC (A145)
- K0223-3 Anti-CXCR7 (RDC1) (Human) mAb (9C4)
- K0223-5 Anti-CXCR7 (RDC1) (Human) mAb-PE (9C4)