

MONOCLONAL ANTIBODY

Anti-HLA-A24 (Human) mAb-PE

Code No.	Clone	Subclass	Quantity
K0208-5	17A10	Mouse IgG2b	1 mL (50 tests)

BACKGROUND: HLA (human leukocyte antigen)-A24 is a class I MHC antigen. HLA-A24 is the most frequent HLA class I molecule in Asian populations, present in approximately ~70% of the Japanese population. HLA-A24 is also found in approximately 35% of the Indian population and 19% of Caucasians. HLA antigens may play a role in genetic susceptibility to disease.

SOURCE: This antibody was purified from hybridoma (clone 17A10) supernatant using protein A agarose. This hybridoma was established by fusion of mouse myeloma cell Sp2/0 with C57BL/6 Tg mouse splenocyte immunized with the recombinant human HLA-A24.

FORMULATION: 50 tests in 1 mL volume of PBS containing 1% BSA and 0.09% NaN₃.

*Azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush plenty of water when disposing materials containing azide into drain.

STORAGE: This antibody solution is stable for one year from the date of purchase when stored at 4°C.

REACTIVITY: This antibody reacts with HLA-A24 on Flow cytometry.

Note: It was reported that this clone 17A10 cross-reacted to HLA-B27 and some indeterminate HLA. Although HLA-B27 population is so small in Japanese, about 20% of tested population in our laboratories reacted to this antibody as false-positive. To ensure your experiment, you should confirm HLA genotyping.

APPLICATION:

Flow cytometry: 20 µL (Ready for use)

*Please refer to the data sheet (MBL; code no. K0208-3) for other applications.

Detailed procedure is provided in the following **PROTOCOLS**.

SPECIES CROSS REACTIVITY:

Species	Human	Mouse	Rat
Cell	LCL721	Not tested	Not tested
Reactivity on FCM	+		


INTENDED USE:

For Research Use Only. Not for use in diagnostic procedures.

REFERENCES:

- 1) Najima, Y., *et al.*, *Blood* **127**, 722-734 (2016) [FCM]
- 2) Kozako, T., *et al.*, *J. Immunol.* **177**, 5718-5726 (2006)
- 3) Lutz, C. T., *et al.*, *J. Immunol.* **153**, 4099-4110 (1994)
- 4) Tahara, T., *et al.*, *Immunogenetics* **32**, 351-360 (1990)

Clone 17A10 is used in these references.



Flow cytometric analysis of HLA-A24 expression on LCL721 cells (right) and Jurkat cells (left). Open histogram indicate the reaction of isotypic control to the cells. Shaded histograms indicate the reaction of K0208-5 to the cells.

PROTOCOLS:

Flow cytometric analysis for floating cells

We usually use Fisher tubes or equivalents as reaction tubes for all steps described below.

- 1) Wash the cells 3 times with washing buffer [PBS containing 2% fetal calf serum (FCS) and 0.09% NaN₃].
*Azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush plenty of water when disposing materials containing azide into drain.
- 2) Resuspend the cells with washing buffer (5x10⁶ cells/mL).
- 3) Add 50 µL of the cell suspension into each tube, and centrifuge at 500 x g for 1 minute at room temperature (20~25°C). Remove supernatant by careful aspiration.
- 4) Add 20 µL of Clear Back (human Fc receptor blocking reagent, MBL; code no. MTG-001) to the cell pellet after tapping. Mix well and incubate for 5 minutes at room temperature.
- 5) Add 20 µL of the primary antibody. Mix well and incubate for 30 minutes at room temperature.

- 6) Add 1 mL of the washing buffer followed by centrifugation at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 7) Resuspend the cells with 500 µL of the washing buffer and analyze by a flow cytometer.

(Positive control for Flow cytometry; LCL721)

Flow cytometric analysis for whole blood cells

We usually use Fisher tubes or equivalents as reaction tubes for all step described below.

- 1) Add 20 µL of the primary antibody into each tube.
- 2) Add 50 µL of whole blood into each tube. Mix well, and incubate for 30 minutes at room temperature (20~25°C).
- 3) Add 1 mL of washing buffer [PBS containing 2% fetal calf serum (FCS) and 0.09% NaN₃] followed by centrifugation at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 4) Lyse with OptiLyse C (for analysis on Beckman Coulter instruments) or OptiLyse B (for analysis on BD instruments), using the procedure recommended in the respective package inserts.
- 5) Add 1 mL of H₂O to each tube and incubate for 10 minutes at room temperature.
- 6) Centrifuge at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 7) Add 1 mL of washing buffer followed by centrifugation at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 8) Resuspend the cells with 500 µL of the washing buffer and analyze by a flow cytometer.

RELATED PRODUCTS:

K0186-3	Anti-HLA-A2 (Human) mAb (BB7.2)
K0186-4	Anti-HLA-A2 (Human) mAb-FITC (BB7.2)
K0186-5	Anti-HLA-A2 (Human) mAb-PE (BB7.2)
K0208-3	Anti-HLA-A24 (Human) mAb (17A10)
K0208-4	Anti-HLA-A24 (Human) mAb-FITC (17A10)
K0208-5	Anti-HLA-A24 (Human) mAb-PE (17A10)
K0208-A48	Anti-HLA-A24 (Human) mAb -Alexa Fluor [®] 488 (17A10)
K0208-A64	Anti-HLA-A24 (Human) mAb -Alexa Fluor [®] 647 (17A10)
K0209-3	Anti-HLA-A24 (Human) mAb (22E1)
K0209-4	Anti-HLA-A24 (Human) mAb-FITC (22E1)
K0209-5	Anti-HLA-A24 (Human) mAb-PE (22E1)
D226-3	Anti-HLA class I (HLA-A,B,C) (Human) mAb (EMR8-5)
K0126-3	Anti-HLA-E (Human) mAb (MEM-E/02)
K0215-3	Anti-HLA-E (Human) mAb (4D12)
K0125-3	Anti-HLA-G (Human) mAb (MEM-G/1)
K0019-1	Anti-HLA-DR (Human) mAb (LN-3)
M077-5	Mouse IgG2b (isotype control)-PE (3D12)