

 **My select** sampler set

Rat IgG1 (isotype control)-Alexa Fluor[®] 488

Code No.	Clone	Subclass	Quantity	Concentration
M080-A48MS	1H5	Rat IgG1	20 μ L	1 mg/mL

SOURCE: This antibody was purified from hybridoma (clone 1H5) supernatant using protein G agarose. This hybridoma was established by fusion of mouse myeloma cell P3U1 with rat lymph nodes immunized with KLH.

FORMULATION: 20 μ g IgG in 20 μ L volume of PBS containing 1% BSA and 0.09% NaN₃.

*Azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush plenty of water when disposing materials containing azide into drain.

STORAGE: This antibody solution is stable for one year from the date of purchase when stored at 4°C.

REACTIVITY: No specific binding is detected on mouse splenocytes.


APPLICATION:

Flow cytometry: This antibody can be used as a negative isotypic control. The concentration will depend on the conditions.

Detailed procedure is provided in the following **PROTOCOL**.

INTENDED USE:

For Research Use Only. Not for use in diagnostic procedures.



Flow cytometric analysis of rat IgG1 reactivity on mouse splenocytes.

PROTOCOL:

Flow cytometric analysis for floating cells

We usually use Fisher tubes or equivalents as reaction tubes for all steps described below.

- 1) Wash the cells 3 times with washing buffer [PBS containing 2% fetal calf serum (FCS) and 0.1% NaN₃].
- 2) Resuspend the cells with washing buffer (6x10⁶ cells/mL).

RELATED PRODUCTS:

Functional grade antibodies

- M080-3M2 Rat IgG1 (isotype control) (1H5)
- M081-3M2 Rat IgG2a (isotype control) (2H3)
- M090-3M2 Rat IgG2b (isotype control) (3G8)
- M075-3M2 Mouse IgG1 (isotype control) (2E12)
- M076-3M2 Mouse IgG2a (isotype control) (6H3)
- M077-3M2 Mouse IgG2b (isotype control) (3D12)

Purified antibodies

- M080-3 Rat IgG1 (isotype control) (1H5)
- M080-4 Rat IgG1 (isotype control)-FITC (1H5)
- M080-5 Rat IgG1 (isotype control)-PE (1H5)
- M080-A48 Rat IgG1 (isotype control)-Alexa Fluor[®] 488 (1H5)
- M081-3 Rat IgG2a (isotype control) (2H3)
- M081-4 Rat IgG2a (isotype control)-FITC (2H3)
- M081-5 Rat IgG2a (isotype control)-PE (2H3)
- M081-8 Rat IgG2a (isotype control)-Agarose (2H3)
- M081-9 Rat IgG2a (isotype control)-Magnetic beads (2H3)
- M081-A48 Rat IgG2a (isotype control)-Alexa Fluor[®] 488 (2H3)
- M090-3 Rat IgG2b (isotype control) (3G8)
- M090-4 Rat IgG2b (isotype control)-FITC (3G8)
- M090-5 Rat IgG2b (isotype control)-PE (3G8)
- M090-A48 Rat IgG2b (isotype control)-Alexa Fluor[®] 488 (3G8)
- M082-3 Rat IgG2c (isotype control) (6E12)
- M082-4 Rat IgG2c (isotype control)-FITC (6E12)
- M075-3 Mouse IgG1 (isotype control) (2E12)
- M075-4 Mouse IgG1 (isotype control)-FITC (2E12)
- M075-5 Mouse IgG1 (isotype control)-PE (2E12)
- M075-8 Mouse IgG1 (isotype control)-Agarose (2E12)
- M075-9 Mouse IgG1 (isotype control)-Magnetic beads (2E12)
- M075-A48 Mouse IgG1 (isotype control)-Alexa Fluor[®] 488 (2E12)
- M075-A64 Mouse IgG1 (isotype control)-Alexa Fluor[®] 647 (2E12)
- M076-3 Mouse IgG2a (isotype control) (6H3)
- M076-4 Mouse IgG2a (isotype control)-FITC (6H3)
- M076-5 Mouse IgG2a (isotype control)-PE (6H3)
- M076-9 Mouse IgG2a (isotype control)-Magnetic beads (6H3)
- M076-A48 Mouse IgG2a (isotype control)-Alexa Fluor[®] 488 (6H3)
- M076-A64 Mouse IgG2a (isotype control)-Alexa Fluor[®] 647 (6H3)
- M077-3 Mouse IgG2b (isotype control) (3D12)
- M077-4 Mouse IgG2b (isotype control)-FITC (3D12)
- M077-5 Mouse IgG2b isotype control-PE (3D12)
- M077-9 Mouse IgG2b (isotype control)-Magnetic beads (3D12)
- M077-A48 Mouse IgG2b (isotype control)-Alexa Fluor[®] 488 (3D12)
- M077-A64 Mouse IgG2b (isotype control)-Alexa Fluor[®] 647 (3D12)
- M078-3 Mouse IgG3 (isotype control) (6A3)
- M078-4 Mouse IgG3 (isotype control)-FITC (6A3)
- M079-3 Mouse IgM (isotype control) (7E10)
- M194-3 Human IgG1 isotype control chimeric mAb
- M195-3 Human IgG2 isotype control chimeric mAb
- PM035 Normal Rabbit IgG (polyclonal)
- PM035-8 Normal Rabbit IgG-Agarose (polyclonal)
- PM067 Normal Guinea Pig IgG (polyclonal)
- M189-3 Syrian Hamster IgG (isotype control)
- M199-3 Armenian Hamster IgG (isotype control)

Other related antibodies and kits are also available.
Please visit our web site at <http://ruo.mbl.co.jp>

130603-1

LABEL LICENSES:

This product is provided under an agreement between LIFE TECHNOLOGIES Corporation, and Medical & Biological Laboratories Co., LTD, and the manufacture, use, sale or import of this product may be subject to one or more U.S. patents, pending applications, and corresponding non-U.S. equivalents, owned by or licensed to LIFE TECHNOLOGIES Corporation. For information on purchasing a license to this product for any other use, contact LIFE TECHNOLOGIES Corporation, Molecular Probes Detection Technologies, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, USA, Tel: (541) 465-8300. Fax: (541) 335-0504.

Alexa Fluor[®] is a registered trademark of Molecular Probes, Inc.