

**For Research Use Only.
Not for use in diagnostic procedures.**



Rat IgG2c isotype control (Functional Grade)

CODE No.	M082-3M2
CLONALITY	Monoclonal
CLONE	6E12
ISOTYPE	Rat IgG2c
QUANTITY	100 µL, 5 mg/mL
SOURCE	50% saturated ammonium sulfate precipitated fraction from mouse ascites fluid
IMMUNOGEN	KLH
FORMURATION	5 mg/mL in PBS. Azide free, 0.22 µm sterile-filtered Endotoxin level is <0.5 EU/mg antibody, as determined by the LAL assay
STORAGE	This antibody solution is stable for one year from the date of purchase when stored at -20°C. Avoid repeated freezing and thawing.

For more information, please visit our web site <https://ruo.mbl.co.jp/>

RELATED PRODUCTS

[Functional grade antibody]

M075-3M2 Mouse IgG1 isotype control FG (2E12)
M076-3M2 Mouse IgG2a isotype control FG (6H3)
M077-3M2 Mouse IgG2b isotype control FG (3D12)
M078-3M2 Mouse IgG3 isotype control FG (6A3)
M080-3M2 Rat IgG1 isotype control FG (1H5)
M081-3M2 Rat IgG2a isotype control FG (2H3)
M090-3M2 Rat IgG2b isotype control FG (3G8)

[Purified antibody]

M075-3 Mouse IgG1 isotype control (2E12)
M075-4 Mouse IgG1 isotype control-FITC (2E12)
M075-5 Mouse IgG1 isotype control-PE (2E12)
M075-8 Mouse IgG1 isotype control-Agarose (2E12)
M075-A48 Mouse IgG1 isotype control-Alexa Fluor[®] 488 (2E12)
M075-A64 Mouse IgG1 isotype control-Alexa Fluor[®] 647 (2E12)
M076-3 Mouse IgG2a isotype control (6H3)
M076-4 Mouse IgG2a isotype control-FITC (6H3)
M076-5 Mouse IgG2a isotype control-PE (6H3)
M076-A48 Mouse IgG2a isotype control-Alexa Fluor[®] 488 (6H3)
M076-A64 Mouse IgG2a isotype control-Alexa Fluor[®] 647 (6H3)
M077-3 Mouse IgG2b isotype control (3D12)
M077-4 Mouse IgG2b isotype control-FITC (3D12)
M077-5 Mouse IgG2b isotype control-PE (3D12)
M077-A48 Mouse IgG2b isotype control-Alexa Fluor[®] 488 (3D12)
M077-A64 Mouse IgG2b isotype control-Alexa Fluor[®] 647 (3D12)
M078-3 Mouse IgG3 isotype control (6A3)
M078-4 Mouse IgG3 isotype control-FITC (6A3)
M079-3 Mouse IgM isotype control (7E10)
M080-3 Rat IgG1 isotype control (1H5)
M080-4 Rat IgG1 isotype control-FITC (1H5)
M080-5 Rat IgG1 isotype control-PE (1H5)
M080-A48 Rat IgG1 isotype control-Alexa Fluor[®] 488 (1H5)
M081-3 Rat IgG2a isotype control (2H3)
M081-4 Rat IgG2a isotype control-FITC (2H3)
M081-5 Rat IgG2a isotype control-PE (2H3)
M081-8 Rat IgG2a isotype control-Agarose (2H3)
M081-A48 Rat IgG2a isotype control-Alexa Fluor[®] 488 (2H3)
M090-3 Rat IgG2b isotype control (3G8)
M090-4 Rat IgG2b isotype control-FITC (3G8)
M090-5 Rat IgG2b isotype control-PE (3G8)
M090-A48 Rat IgG2b isotype control-Alexa Fluor[®] 488 (3G8)
M082-3 Rat IgG2c isotype control (6E12)
M082-4 Rat IgG2c isotype control-FITC (6E12)
PM035 Normal Rabbit IgG (polyclonal)
PM035-8 Normal Rabbit IgG-Agarose (polyclonal)