

Announcement: New Arrivals

New Arrivals of MBL products, 28 items

Dear Customers,

First of all, please accept our sincere gratitude for your continued support.

We are pleased to inform you that we are launching small size products for particularly popular items from the anti-epitope tag antibody series. It is very useful when you want to do first time trial, compare with several antibodies, or need only a small volume. Please try them this opportunity.

Sincerely,
MBL
May, 2020

[Effective date]

June 1st, 2020

[Product information]

Please check Appendix

<Appendix>

Code No.	Product	Size
561MS	Anti-HA-tag pAb for My select sampler set	20 µL
562MS	Anti-Myc-tag pAb for My select sampler set	20 µL
598MS	Anti-GFP pAb for My select sampler set	20 µL
D153-11MS	Anti-GFP mAb-Magnetic Beads for My select sampler set	4 tests (Slurry: 200 µL)
D153-3MS	Anti-GFP mAb for My select sampler set	20 µg/20 µL
D291-11MS	Anti-His-tag mAb-Magnetic Beads for My select sampler set	4 tests (Slurry: 200 µL)
D291-3MS	Anti-His-tag mAb for My select sampler set	20 µg/20 µL
M047-11MS	Anti-Myc-tag mAb-Magnetic Beads for My select sampler set	4 tests (Slurry: 200 µL)
M047-3MS	Anti-Myc-tag mAb for My select sampler set	40 µg/20 µL
M048-3MS	Anti-GFP mAb for My select sampler set	20 µg/20 µL
M089-3MS	Anti-His-tag mAb for My select sampler set	20 µg/20 µL
M132-3MS	Anti-HA-tag mAb for My select sampler set	40 µg/20 µL
M136-3MS	Anti-His-tag mAb for My select sampler set	20 µg/20 µL
M155-3MS	Anti-RFP mAb for My select sampler set	20 µg/20 µL
M165-3MS	Anti-RFP mAb for My select sampler set	20 µg/20 µL
M167-3MS	Anti-V5-tag mAb for My select sampler set	20 µg/20 µL
M180-11MS	Anti-HA-tag mAb-Magnetic Beads for My select sampler set	4 tests (Slurry: 200 µL)
M180-3MS	Anti-HA-tag mAb for My select sampler set	20 µg/20 µL
M185-11MS	Anti-DDDDK-tag mAb-Magnetic Beads for My select sampler set	4 tests (Slurry: 200 µL)
M185-3MS	Anti-DDDDK-tag mAb for My select sampler set	20 µg/20 µL
M192-3MS	Anti-Myc-tag mAb for My select sampler set	20 µg/20 µL
M204-3MS	Anti-RFP mAb for My select sampler set	20 µg/20 µL
M215-11MS	Anti-V5-tag mAb-Magnetic Beads for My select sampler set	4 tests (Slurry: 200 µL)
M215-3MS	Anti-V5-tag mAb for My select sampler set	20 µg/20 µL
PM003MS	Anti-V5-tag pAb for My select sampler set	20 µL
PM005MS	Anti-RFP pAb for My select sampler set	20 µL
PM020MS	Anti-DDDDK-tag pAb for My select sampler set	20 µL
PM032MS	Anti-His-tag pAb for My select sampler set	20 µL